

< e-CHECK(XS) ASSAY SHEET >

For XS-Series

Lot No. 61730804

Exp. date 2026-09-13

L1 : Level 1

	RBC 10 ⁶ /uL	HGB		HCT %	MCV fL	MCH		MCHC		PLT 10 ³ /uL
		g/dL	mmol/L			pg	fmol	g/dL	mmol/L	
TARGET	2.56	5.6	3.5	18.0	70.5	21.9	1.36	31.1	19.3	60
LIMIT	0.13	0.2	0.1	1.1	3.5	1.1	0.07	1.9	1.2	24

	RDW-SD fL	RDW-CV %	PDW fL	PCT %	MPV fL	P-LCR %
TARGET	46.2	18.4	9.0	0.06	9.6	17.0
LIMIT	6.9	2.8	3.6	0.03	1.4	14.5

	WBC-D 10 ³ /uL	WBC-C 10 ³ /uL	NEUT# 10 ³ /uL	LYMPH# 10 ³ /uL	MONO# 10 ³ /uL	EO# 10 ³ /uL	BASO# 10 ³ /uL
TARGET	2.87	2.77	1.10	1.07	0.27	0.27	0.16
LIMIT	0.29	0.28	0.33	0.43	0.24	0.14	0.12

Please refer to the help message provided in the assay CD-ROM for the handling of the e-CHECK(XS).

Please make sure the lot number before reading the data from the assay CD-ROM.

Do not leave control blood in the room temperature over an hour.

Please store it in a refrigerator (2-8°C) immediately after use.

U.S. Patents 5,432,089; 6,200,500; 6,221,668; 6,399,388; 6,403,377; 6,406,915

	DIFF-X	DIFF-Y	FSC-X
TARGET	169.9	66.9	27.0
LIMIT	15.0	30.0	15.0

**Confidential
Service Staff Only**

< e-CHECK(XS) ASSAY SHEET >

For XS-Series

Lot No. 61730805

Exp. date 2026-09-13

L2 : Level 2

	RBC	HGB		HCT	MCV	MCH		MCHC		PLT
	10 ⁶ /uL	g/dL	mmol/L	%	fL	pg	fmol	g/dL	mmol/L	10 ³ /uL
TARGET	4.24	12.1	7.5	36.3	85.6	28.5	1.77	33.3	20.7	220
LIMIT	0.21	0.4	0.2	1.8	4.3	1.4	0.09	2.0	1.2	33

	RDW-SD	RDW-CV	PDW	PCT	MPV	P-LCR
	fL	%	fL	%	fL	%
TARGET	49.8	16.2	9.2	0.22	10.0	17.6
LIMIT	7.5	2.4	2.8	0.07	1.5	8.8

	WBC-D	WBC-C	NEUT#	LYMPH#	MONO#	EO#	BASO#
	10 ³ /uL	10 ³ /uL	10 ³ /uL	10 ³ /uL	10 ³ /uL	10 ³ /uL	10 ³ /uL
TARGET	6.73	6.49	2.80	2.16	0.69	0.67	0.41
LIMIT	0.40	0.39	0.84	0.65	0.62	0.34	0.32

Please refer to the help message provided in the assay CD-ROM for the handling of the e-CHECK(XS).

Please make sure the lot number before reading the data from the assay CD-ROM.

Do not leave control blood in the room temperature over an hour.

Please store it in a refrigerator (2-8°C) immediately after use.

U.S. Patents 5,432,089; 6,200,500; 6,221,668; 6,399,388; 6,403,377; 6,406,915

	DIFF-X	DIFF-Y	FSC-X
TARGET	169.9	69.1	27.0
LIMIT	15.0	30.0	15.0

**Confidential
Service Staff Only**

< e-CHECK(XS) ASSAY SHEET >

For XS-Series

Lot No. 61730806

Exp. date 2026-09-13

L3 : Level 3

	RBC 10 ⁶ /uL	HGB		HCT %	MCV fL	MCH		MCHC		PLT 10 ³ /uL
		g/dL	mmol/L			pg	fmol	g/dL	mmol/L	
TARGET	5.10	15.8	9.8	47.1	92.3	31.0	1.92	33.5	20.8	524
LIMIT	0.26	0.5	0.3	2.4	4.6	1.6	0.10	2.0	1.2	68

	RDW-SD fL	RDW-CV %	PDW fL	PCT %	MPV fL	P-LCR %
TARGET	47.2	14.9	9.7	0.54	10.3	19.7
LIMIT	7.1	2.2	1.9	0.14	1.5	9.9

	WBC-D 10 ³ /uL	WBC-C 10 ³ /uL	NEUT# 10 ³ /uL	LYMPH# 10 ³ /uL	MONO# 10 ³ /uL	EO# 10 ³ /uL	BASO# 10 ³ /uL
TARGET	16.09	15.50	7.20	4.31	1.73	1.74	1.11
LIMIT	0.97	0.93	2.16	1.29	1.56	0.87	0.87

Please refer to the help message provided in the assay CD-ROM for the handling of the e-CHECK(XS).

Please make sure the lot number before reading the data from the assay CD-ROM.

Do not leave control blood in the room temperature over an hour.

Please store it in a refrigerator (2-8°C) immediately after use.

U.S. Patents 5,432,089; 6,200,500; 6,221,668; 6,399,388; 6,403,377; 6,406,915

	DIFF-X	DIFF-Y	FSC-X
TARGET	169.2	60.7	27.6
LIMIT	15.0	30.0	15.0

**Confidential
Service Staff Only**