

< e-CHECK(XS) ASSAY SHEET >

For XS-Series

Lot No. 60610804

Exp. date 2026-05-24

L1 : Level 1

	RBC 10 ⁶ /uL	HGB		HCT %	MCV fL	MCH		MCHC		PLT 10 ³ /uL
		g/dL	mmol/L			pg	fmol	g/dL	mmol/L	
TARGET	2.60	6.3	3.9	19.5	74.9	24.2	1.50	32.3	20.1	59
LIMIT	0.13	0.3	0.2	1.2	3.7	1.2	0.08	1.9	1.2	24

	RDW-SD fL	RDW-CV %	PDW fL	PCT %	MPV fL	P-LCR %
TARGET	51.2	19.1	8.4	0.06	9.3	13.8
LIMIT	7.7	2.9	3.4	0.03	1.4	11.7

	WBC-D 10 ³ /uL	WBC-C 10 ³ /uL	NEUT# 10 ³ /uL	LYMPH# 10 ³ /uL	MONO# 10 ³ /uL	EO# 10 ³ /uL	BASO# 10 ³ /uL
TARGET	2.94	2.84	1.15	1.06	0.29	0.27	0.17
LIMIT	0.29	0.28	0.35	0.42	0.26	0.14	0.13

Please refer to the help message provided in the assay CD-ROM for the handling of the e-CHECK(XS).

Please make sure the lot number before reading the data from the assay CD-ROM.

Do not leave control blood in the room temperature over an hour.

Please store it in a refrigerator (2-8°C) immediately after use.

U.S. Patents 5,432,089; 6,200,500; 6,221,668; 6,399,388; 6,403,377; 6,406,915

	DIFF-X	DIFF-Y	FSC-X
TARGET	169.6	67.0	26.7
LIMIT	15.0	30.0	15.0

Confidential
Service Staff Only

< e-CHECK(XS) ASSAY SHEET >

For XS-Series

Lot No. 60610805

Exp. date 2026-05-24

L2 : Level 2

	RBC 10 ⁶ /uL	HGB		HCT %	MCV fL	MCH		MCHC		PLT 10 ³ /uL
		g/dL	mmol/L			pg	fmol	g/dL	mmol/L	
TARGET	4.31	12.0	7.4	36.6	85.0	27.8	1.73	32.8	20.4	220
LIMIT	0.22	0.4	0.2	1.8	4.3	1.4	0.09	2.0	1.2	33

	RDW-SD fL	RDW-CV %	PDW fL	PCT %	MPV fL	P-LCR %
TARGET	47.4	15.7	8.6	0.21	9.7	14.7
LIMIT	7.1	2.4	2.6	0.06	1.5	7.4

	WBC-D 10 ³ /uL	WBC-C 10 ³ /uL	NEUT# 10 ³ /uL	LYMPH# 10 ³ /uL	MONO# 10 ³ /uL	EO# 10 ³ /uL	BASO# 10 ³ /uL
TARGET	6.74	6.53	2.90	2.09	0.66	0.66	0.43
LIMIT	0.40	0.39	0.87	0.63	0.59	0.33	0.34

Please refer to the help message provided in the assay CD-ROM for the handling of the e-CHECK(XS).

Please make sure the lot number before reading the data from the assay CD-ROM.

Do not leave control blood in the room temperature over an hour.

Please store it in a refrigerator (2-8°C) immediately after use.

U.S. Patents 5,432,089; 6,200,500; 6,221,668; 6,399,388; 6,403,377; 6,406,915

	DIFF-X	DIFF-Y	FSC-X
TARGET	166.8	60.9	28.3
LIMIT	15.0	30.0	15.0

Confidential
Service Staff Only

< e-CHECK(XS) ASSAY SHEET >

For XS-Series

Lot No. 60610806

Exp. date 2026-05-24

L3 : Level 3

	RBC 10 ⁶ /uL	HGB		HCT %	MCV fL	MCH		MCHC		PLT 10 ³ /uL
		g/dL	mmol/L			pg	fmol	g/dL	mmol/L	
TARGET	5.20	15.8	9.8	47.1	90.6	30.4	1.89	33.5	20.8	522
LIMIT	0.26	0.5	0.3	2.4	4.5	1.5	0.09	2.0	1.2	68

	RDW-SD fL	RDW-CV %	PDW fL	PCT %	MPV fL	P-LCR %
TARGET	49.9	16.0	9.0	0.52	9.9	16.3
LIMIT	7.5	2.4	1.8	0.13	1.5	8.2

	WBC-D 10 ³ /uL	WBC-C 10 ³ /uL	NEUT# 10 ³ /uL	LYMPH# 10 ³ /uL	MONO# 10 ³ /uL	EO# 10 ³ /uL	BASO# 10 ³ /uL
TARGET	16.52	16.02	7.67	4.20	1.72	1.77	1.16
LIMIT	0.99	0.96	2.30	1.26	1.55	0.89	0.90

Please refer to the help message provided in the assay CD-ROM for the handling of the e-CHECK(XS).

Please make sure the lot number before reading the data from the assay CD-ROM.

Do not leave control blood in the room temperature over an hour.

Please store it in a refrigerator (2-8°C) immediately after use.

U.S. Patents 5,432,089; 6,200,500; 6,221,668; 6,399,388; 6,403,377; 6,406,915

	DIFF-X	DIFF-Y	FSC-X
TARGET	166.9	61.1	28.7
LIMIT	15.0	30.0	15.0

Confidential
Service Staff Only