

< e-CHECK(XS) ASSAY SHEET >

For XS-Series

Lot No. 62290804

Exp. date 2026-11-08

L1 : Level 1

	RBC 10 ⁶ /uL	HGB		HCT %	MCV fL	MCH		MCHC		PLT 10 ³ /uL
		g/dL	mmol/L			pg	fmol	g/dL	mmol/L	
TARGET	2.56	6.0	3.7	18.7	73.1	23.4	1.45	32.1	19.9	58
LIMIT	0.13	0.2	0.1	1.1	3.7	1.2	0.07	1.9	1.2	23

	RDW-SD fL	RDW-CV %	PDW fL	PCT %	MPV fL	P-LCR %
TARGET	46.3	17.4	9.8	0.06	9.8	17.9
LIMIT	6.9	2.6	3.9	0.03	1.5	15.2

	WBC-D 10 ³ /uL	WBC-C 10 ³ /uL	NEUT# 10 ³ /uL	LYMPH# 10 ³ /uL	MONO# 10 ³ /uL	EO# 10 ³ /uL	BASO# 10 ³ /uL
TARGET	3.05	2.94	1.13	1.17	0.30	0.28	0.17
LIMIT	0.31	0.29	0.34	0.47	0.27	0.14	0.13

Please refer to the help message provided in the assay CD-ROM for the handling of the e-CHECK(XS).

Please make sure the lot number before reading the data from the assay CD-ROM.

Do not leave control blood in the room temperature over an hour.

Please store it in a refrigerator (2-8°C) immediately after use.

U.S. Patents 5,432,089; 6,200,500; 6,221,668; 6,399,388; 6,403,377; 6,406,915

	DIFF-X	DIFF-Y	FSC-X
TARGET	168.8	64.3	27.1
LIMIT	15.0	30.0	15.0

**Confidential
Service Staff Only**

< e-CHECK(XS) ASSAY SHEET >

For XS-Series

Lot No. 62290805

Exp. date 2026-11-08

L2 : Level 2

	RBC 10 ⁶ /uL	HGB		HCT %	MCV fL	MCH		MCHC		PLT 10 ³ /uL
		g/dL	mmol/L			pg	fmol	g/dL	mmol/L	
TARGET	4.25	12.1	7.5	36.3	85.3	28.5	1.77	33.3	20.7	214
LIMIT	0.21	0.4	0.2	1.8	4.3	1.4	0.09	2.0	1.2	32

	RDW-SD fL	RDW-CV %	PDW fL	PCT %	MPV fL	P-LCR %
TARGET	45.2	14.9	9.8	0.22	10.0	18.3
LIMIT	6.8	2.2	2.9	0.07	1.5	9.2

	WBC-D 10 ³ /uL	WBC-C 10 ³ /uL	NEUT# 10 ³ /uL	LYMPH# 10 ³ /uL	MONO# 10 ³ /uL	EO# 10 ³ /uL	BASO# 10 ³ /uL
TARGET	7.13	6.85	2.86	2.41	0.74	0.68	0.44
LIMIT	0.43	0.41	0.86	0.72	0.67	0.34	0.34

Please refer to the help message provided in the assay CD-ROM for the handling of the e-CHECK(XS).

Please make sure the lot number before reading the data from the assay CD-ROM.

Do not leave control blood in the room temperature over an hour.

Please store it in a refrigerator (2-8°C) immediately after use.

U.S. Patents 5,432,089; 6,200,500; 6,221,668; 6,399,388; 6,403,377; 6,406,915

	DIFF-X	DIFF-Y	FSC-X
TARGET	167.7	60.0	26.6
LIMIT	15.0	30.0	15.0

**Confidential
Service Staff Only**

< e-CHECK(XS) ASSAY SHEET >

For XS-Series

Lot No. 62290806

Exp. date 2026-11-08

L3 : Level 3

	RBC 10 ⁶ /uL	HGB		HCT %	MCV fL	MCH		MCHC		PLT 10 ³ /uL
		g/dL	mmol/L			pg	fmol	g/dL	mmol/L	
TARGET	5.12	15.7	9.7	46.2	90.2	30.7	1.90	34.0	21.1	504
LIMIT	0.26	0.5	0.3	2.3	4.5	1.5	0.10	2.0	1.3	66

	RDW-SD fL	RDW-CV %	PDW fL	PCT %	MPV fL	P-LCR %
TARGET	43.5	13.9	10.1	0.52	10.3	20.0
LIMIT	6.5	2.1	2.0	0.13	1.5	10.0

	WBC-D 10 ³ /uL	WBC-C 10 ³ /uL	NEUT# 10 ³ /uL	LYMPH# 10 ³ /uL	MONO# 10 ³ /uL	EO# 10 ³ /uL	BASO# 10 ³ /uL
TARGET	16.25	15.63	7.15	4.57	1.73	1.70	1.10
LIMIT	0.98	0.94	2.15	1.37	1.56	0.85	0.86

Please refer to the help message provided in the assay CD-ROM for the handling of the e-CHECK(XS).

Please make sure the lot number before reading the data from the assay CD-ROM.

Do not leave control blood in the room temperature over an hour.

Please store it in a refrigerator (2-8°C) immediately after use.

U.S. Patents 5,432,089; 6,200,500; 6,221,668; 6,399,388; 6,403,377; 6,406,915

	DIFF-X	DIFF-Y	FSC-X
TARGET	169.9	62.4	26.9
LIMIT	15.0	30.0	15.0

**Confidential
Service Staff Only**