

< e-CHECK(XS) ASSAY SHEET >

For XS-Series

Lot No. 61170804

Exp. date 2026-07-19

L1 : Level 1

	RBC 10 ⁶ /uL	HGB		HCT %	MCV fL	MCH		MCHC		PLT 10 ³ /uL
		g/dL	mmol/L			pg	fmol	g/dL	mmol/L	
TARGET	2.56	6.1	3.8	19.1	74.8	23.8	1.48	31.9	19.8	59
LIMIT	0.13	0.2	0.2	1.1	3.7	1.2	0.07	1.9	1.2	24

	RDW-SD fL	RDW-CV %	PDW fL	PCT %	MPV fL	P-LCR %
TARGET	46.9	17.3	9.0	0.06	10.0	19.1
LIMIT	7.0	2.6	3.6	0.03	1.5	16.2

	WBC-D 10 ³ /uL	WBC-C 10 ³ /uL	NEUT# 10 ³ /uL	LYMPH# 10 ³ /uL	MONO# 10 ³ /uL	EO# 10 ³ /uL	BASO# 10 ³ /uL
TARGET	3.03	2.93	1.17	1.10	0.29	0.29	0.18
LIMIT	0.30	0.29	0.35	0.44	0.26	0.15	0.14

Please refer to the help message provided in the assay CD-ROM for the handling of the e-CHECK(XS).

Please make sure the lot number before reading the data from the assay CD-ROM.

Do not leave control blood in the room temperature over an hour.

Please store it in a refrigerator (2-8°C) immediately after use.

U.S. Patents 5,432,089; 6,200,500; 6,221,668; 6,399,388; 6,403,377; 6,406,915

	DIFF-X	DIFF-Y	FSC-X
TARGET	167.5	69.7	27.9
LIMIT	15.0	30.0	15.0

Confidential
Service Staff Only

< e-CHECK(XS) ASSAY SHEET >

For XS-Series

Lot No. 61170805

Exp. date 2026-07-19

L2 : Level 2

	RBC 10 ⁶ /uL	HGB		HCT %	MCV fL	MCH		MCHC		PLT 10 ³ /uL
		g/dL	mmol/L			pg	fmol	g/dL	mmol/L	
TARGET	4.26	11.6	7.2	34.9	81.9	27.2	1.69	33.2	20.6	220
LIMIT	0.21	0.3	0.2	1.7	4.1	1.4	0.08	2.0	1.2	33

	RDW-SD fL	RDW-CV %	PDW fL	PCT %	MPV fL	P-LCR %
TARGET	46.1	15.8	9.2	0.23	10.3	19.8
LIMIT	6.9	2.4	2.8	0.07	1.5	9.9

	WBC-D 10 ³ /uL	WBC-C 10 ³ /uL	NEUT# 10 ³ /uL	LYMPH# 10 ³ /uL	MONO# 10 ³ /uL	EO# 10 ³ /uL	BASO# 10 ³ /uL
TARGET	6.70	6.47	2.79	2.12	0.68	0.70	0.41
LIMIT	0.40	0.39	0.84	0.64	0.61	0.35	0.32

Please refer to the help message provided in the assay CD-ROM for the handling of the e-CHECK(XS).

Please make sure the lot number before reading the data from the assay CD-ROM.

Do not leave control blood in the room temperature over an hour.

Please store it in a refrigerator (2-8°C) immediately after use.

U.S. Patents 5,432,089; 6,200,500; 6,221,668; 6,399,388; 6,403,377; 6,406,915

	DIFF-X	DIFF-Y	FSC-X
TARGET	168.9	65.5	27.8
LIMIT	15.0	30.0	15.0

Confidential
Service Staff Only

< e-CHECK(XS) ASSAY SHEET >

For XS-Series

Lot No. 61170806

Exp. date 2026-07-19

L3 : Level 3

	RBC 10 ⁶ /uL	HGB		HCT %	MCV fL	MCH		MCHC		PLT 10 ³ /uL
		g/dL	mmol/L			pg	fmol	g/dL	mmol/L	
TARGET	5.20	16.1	10.0	47.2	90.8	31.0	1.92	34.1	21.2	537
LIMIT	0.26	0.5	0.3	2.4	4.5	1.6	0.10	2.0	1.3	70

	RDW-SD fL	RDW-CV %	PDW fL	PCT %	MPV fL	P-LCR %
TARGET	45.1	14.2	9.8	0.56	10.5	21.6
LIMIT	6.8	2.1	2.0	0.14	1.6	10.8

	WBC-D 10 ³ /uL	WBC-C 10 ³ /uL	NEUT# 10 ³ /uL	LYMPH# 10 ³ /uL	MONO# 10 ³ /uL	EO# 10 ³ /uL	BASO# 10 ³ /uL
TARGET	16.37	15.79	7.38	4.38	1.76	1.72	1.13
LIMIT	0.98	0.95	2.21	1.31	1.58	0.86	0.88

Please refer to the help message provided in the assay CD-ROM for the handling of the e-CHECK(XS).

Please make sure the lot number before reading the data from the assay CD-ROM.

Do not leave control blood in the room temperature over an hour.

Please store it in a refrigerator (2-8°C) immediately after use.

U.S. Patents 5,432,089; 6,200,500; 6,221,668; 6,399,388; 6,403,377; 6,406,915

	DIFF-X	DIFF-Y	FSC-X
TARGET	169.4	64.5	27.4
LIMIT	15.0	30.0	15.0

Confidential
Service Staff Only